

Boston Scientific


Advancing science for life™



VERCISE™ CHARGING SYSTEM

QUICK START GUIDE

1. PLUGGING IN THE POWER SUPPLY




Plug the Power Supply into an outlet and connect to the Base Station.

2. CHARGING AND RE-CHARGING

Charge the Charger by placing it in the Base Station.



When your Charger is fully charged, the light will turn green.



3. CHARGING YOUR STIMULATOR:

You have two charging options.

Option A: Using the Charging Collar



Press the Power Button.

Place the Charger into the appropriate pocket on the Charging Collar with the Power Button facing out.

Place the Counter-weight in the pocket opposite of the Charger.

Put the Charging Collar over your neck such that the Charger is centered over the Stimulator.

You can locate the Stimulator by moving the Charger around until it stops beeping.*

Option B: Using the Adhesive Patch



Peel off the clear plastic liner and affix the patch to your Charger.

Press the Power Button.

Locate the Stimulator by moving the Charger around until it stops beeping.*

Peel off the beige paper liner and affix your charger over the implant site.

4. CHARGING COMPLETE

The Charger will emit a double beep when the Stimulator is fully charged.

Return the Charger to the Base Station.

5. SEE HOW IT'S DONE

Visit the following link to view a video on how to charge your Vercise DBS Stimulator:
<https://youtu.be/iDepGIKjxcI>

*The beeping will stop when the charger is partially aligned with the stimulator. See the Vercise™ DBS System Charging Handbook for instructions on how to ensure proper alignment and optimal charging.

Not intended as a substitute for the provided Vercise™ Deep Brain Stimulation (DBS) System Patient Manual or Charging Handbook. The contents of the handbook must be fully understood before using your Vercise Charging System.

Indications for Use: The Vercise™ Deep Brain Stimulation (DBS) System is indicated for use in bilateral stimulation of the subthalamic nucleus (STN) as an adjunctive therapy in reducing some of the symptoms of moderate to advanced levodopa-responsive Parkinson's disease (PD) that are not adequately controlled with medication.

Contraindications: The Vercise DBS System is not recommended for patients who will be exposed to the following procedures: Diathermy as either a treatment for a medical condition or as part of a surgical procedure, Electroconvulsive Therapy (ECT) and Transcranial Magnetic Stimulation (TMS). The safety of these therapies in patients implanted with the Vercise DBS System has not been established. Patients implanted with the Vercise DBS System should not have Magnetic Resonance Imaging (MRI). The Vercise DBS System is not recommended for patients who are unable to operate the system or are poor surgical candidates or who experience unsuccessful test stimulation.

Warnings: Unauthorized modification to the medical devices is prohibited. You should not be exposed to high stimulation levels. High level of stimulation may damage brain tissue. Patients implanted with the Vercise DBS System may be at risk for intracranial hemorrhages (bleeding in the brain) during DBS lead placement. Strong electromagnetic fields, such as power generators, security screeners or theft detections systems, can potentially turn the stimulator off, or cause unpredictable changes in stimulation. The system should not be charged while sleeping. If you notice changes in mood or behavior or have thoughts of suicide contact your physician immediately. Chemical burns may result if the Vercise Stimulator housing is ruptured or pierced. The Deep Brain Stimulation System may interfere with the operation of implanted stimulation devices, such as cardiac pacemakers, implanted cardioverter defibrillators, or medication delivery pumps. Patients should operate motorized vehicles or potentially dangerous machinery with caution. It is unknown if the device may hurt an unborn baby. Your doctor may be able to provide additional information on the Boston Scientific Vercise DBS System. For complete indications for use, contraindications, warnings, precautions, and side effects, call 833-DBS-INFO or 833-327-4636.

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.

All trademarks are the property of their respective owners.

Boston Scientific

Advancing science for life™

25155 Rye Canyon Loop
Valencia, CA 91355 USA

Copyright ©2017 by
Boston Scientific Corporation
or its affiliates. All rights reserved.

www.bostonscientific.com

Printed in the United States.

NM-511502-AA_DEC2017